


Examination Gown




Instructions For Use:

- Peel off the bag, take out the examination gown.
- Double hand into the sleeve. Draw back from the chest.
- Fasten the neck and the waist belt.
- It is use to cover the body to easy examination operations.

Features:

- Material: polypropylene or Spunbond + Meltblown + Spunbond Nonwovens.
- Size: S, M, L, XL, XXL.
- Long sleeve with elastic cuffs.
- Waist belts with ties neck closures.
- Color: blue, green, yellow, white.

Size	Ref	Material
S	229-S-18A	Spunbond + Meltblown + Spunbond Nonwovens
M	229-M-18A	Spunbond + Meltblown + Spunbond Nonwovens
L	229-L-18A	Spunbond + Meltblown + Spunbond Nonwovens
XL	229-XL-18A	Spunbond + Meltblown + Spunbond Nonwovens
XXL	229-XXL-18A	Spunbond + Meltblown + Spunbond Nonwovens
S	229-S-18B	Polypropylene
M	229-M-18B	Polypropylene
L	229-L-18B	Polypropylene
XL	229-XL-18B	Polypropylene
XXL	229-XXL-18B	Polypropylene
S	229-S-18C	Polypropylene
M	229-M-18C	Polypropylene
L	229-L-18C	Polypropylene
XL	229-XL-18C	Polypropylene
XXL	229-XXL-18C	Polypropylene



Caution:

1. This product is for single use only, discard after single use.
2. Do not use it if the package be affected with damp, or the product went moldy.
3. Do not use after expiry date. If over the expiry date, please discard it.
4. Single use only.

Storage:

Store in a cool, dry and dark place. Avoid direct sunlight.

Examination Gown



Secured Medical Direction UK Co., Ltd.

Wayson Commercial Bldg 28 Con-naught Rd West Sheung Wan HK, P.R.C.



Secured Medical Direction UK Co., Ltd.

St. James House 13 Kensington Square W8 5HD London U.K.

www.smd-medical.com